




# Mycetoma Grains Culture Technique



The Mycetoma Research Center,  
University of Khartoum  
WHO Collaborating Center  
on Mycetoma & Skin NTDs

**Nema Ahmed EL Faki**

**Ahmed Fahal**

**The Mycetoma Research Center**

---

# Introduction


Mycetoma is chronic deep implantation mycosis characterised by

- Large painless subcutaneous swellings
- The formation of sinuses,
- Discharge that contains grains.

The frequently located are foot and hand.




# Mycetoma Classification



# Actinomycetoma

## Causative microorganisms

- *Actinomadura madurae*
- *Actinomadura pelletierii*
- *Nocardia brasiliensis*
- *Streptomyces somaliensis*



*Actinomadura pelletierii*




# Eumycetoma

## Causative microorganisms

The Common Causative Agents:

- *Madurella mycetomatis*
- *Trematosphaeria grisea*
- *Scedosporium spp*
- *Acremonium spp*
- *Fusarium spp*



***Madurella mycetomatis***

# Samples Collection & Preparation

Grains from a deep surgical Biopsy must be collected

Grains from open sinuses discharge should be avoided as they are frequently dead and infected.





## Samples Collection & Preparation



Grains from a deep surgical Biopsy must be collected

## Samples Collection & Preparation

Grains should be placed in a sterile container filled with normal saline and labelled with the patient identification and date of collection.






---

## Samples Collection & Preparation

Grains should be washed with normal saline until they appear clear and without any tissue debris.



---

## Samples Collection & Preparation

Eumycetoma grains usually appear as granules or aggregates with a chalky consistency and various colours ranging from brown to black.



---

## Samples Collection & Preparation

Actinomycetoma grains usually appear as granules or aggregates with a chalky consistency and various colours ranging from yellow to red.



---

## Culture Media

For Actinomycetes the Commonly used media are:

- **Blood agar:**  
To support the growth of most bacterial pathogens, including actinomycetes.
- **Yeast extract Agar:**  
To support the growth of slow-growing actinomycetes.
- **Lowenstein Jensen (LJ) Media**






---

## Culture Media

For Eumycetoma, the Commonly used medium is :

### **Sabouraud Dextrose agar**

supplemented with gentamycin  
to inhibit the growth of bacteria



## Culture Technique

The microbiological sterile loop or syringe is used to inoculate the grains into Media at the center of the universal tube.

Culture is repeated in duplicate.



---

## Culture Technique

In the absence of grains, the biopsy is fragmented, rinsed, and then inoculated into the universal tube as sections.



---

## Culture Technique

Lable the universal tube with patient identification and date of inoculation.





## Culture Technique

Document the specimen identification and date of inoculation in the documentation log.



## Culture Technique

Incubate the cultured grains at 25-37°C for 7-10 days for both actinomycetoma and eumycetoma.



---

## Culture Technique

Daily observation for growth is mandatory.



## Actinomycetoma Macroscopical & Microscopical Examination

Upon positive growth being observed, macroscopic identification of colonial morphology must be done.

Microscopic examination using:

- Gram stain
- Modified Ziehl-Neelsen

. (Gram Stain & ZN Stain SOPs in MRC website )






## Actinomycetoma Macroscopical & Microscopical Examination

Upon positive growth being observed, macroscopic identification of colonial morphology must be done.

Microscopic examination using:

- Gram stain
- Modified Ziehl-Neelsen

. (Gram Stain & ZN Stain SOPs in MRC website )



---

## Eumycetoma Macroscopical & Microscopical Examination

Once fungal colonies have developed, they should be examined macroscopically for characteristic features such as colour, texture, and growth pattern.


### **Microscopical Examination:**

microscopic morphology of hyphae and spores are observed using LPCB Stain



(LPCB Stain SOPs in MRC website )


## Eumycetoma Macroscopical & Microscopical Examination



## Subculturing Technique

### For actinomycetoma

1. Once colonies are observed,
2. Subculture them onto fresh media to obtain pure cultures for further characterisation and identification.
3. Using a sterile loop or syringe, gather multiple pieces of the growth and introduce them into yeast extract Agar plate and slope
4. Wrap the plate well to Avoid agar drying.






## Subculturing Technique

### For actinomycetoma

5. Incubate the tubes and plates at 37°C for 3-7 days.

6. Regularly monitor the plates for the development of characteristic colonies



## Subculturing Technique

### For Eumycetoma

1. Collect multiple pieces from growth using a mycological loop or syringe, inoculate into a Sabouraud dextrose plate, and slop.
2. Wrap the plate well to Avoid agar drying
3. Incubate the tubes at 37°C for 7-10 days.
4. Incubate the tubes and plates at 37°C for 3-7 days.
5. Regularly monitor the plates and tubes for the development of characteristic colonies.



## Culture Preservation

1. Preserve pure cultures of the isolated strains for future reference and research in cryotube filled with 10% glycerol in the appropriate storage conditions, commonly in a -80C refrigerator.
2. For Actinomycetes isolates the cryotube filled with peptone water.
3. Culture and grains are stored at a -80C refrigerator for research and shipment purposes.
4. Eumycetoma culture can be stored in glycerol at room temperature



(preservation and shipment SOPS at MRC website )



The Mycetoma Research Center,  
University of Khartoum  
WHO Collaborating Center  
on Mycetoma & Skin NTDs

[www.mycetoma.edu.sd](http://www.mycetoma.edu.sd)